VisuCal-F Frozen Calibrator Plasma

Product #: FRSP0125
FRSP0105

Store at -60°C or lower.

For *In Vitro* Diagnostic Use.

Intended Use

ENGLISH

The VisuCal-F Frozen Calibrator plasma is intended for use in the calibration of coagulation and fibrinolysis assays.

Summary and Principle

The VisuCal-F Frozen Calibrator plasma is an assayed reference plasma prepared from a pool of citrated human plasma collected from a minimum of 20 donors, buffered with 0.02M HEPES buffer, dispensed and rapidly frozen. This plasma can be used to construct a reference curve when measuring the following parameters: fibrinogen (Clauss Method), coagulation factors II, V, VII, VIII, IX, X, XI and XII, antithrombin activity, α 2-antiplasmin, plasminogen, Protein C activity and antigen, Protein S activity and total antigen and von Willebrand Factor antigen and Ristocetin Cofactor.

Reagents

A. Description:

Calibrator Plasma:

Box containing 25 x 1 mL vials of frozen plasma, Prod. # FRSP0125

OR

Box containing 5 x 1 mL vials of frozen plasma, Prod. # FRSP0105

Each box of VisuCal-F Frozen Calibrator Plasma contains a lot-specific Certificate of Analysis with the assigned reference values for each parameter. The reported reference values were determined against, where available, the SSC/ISTH secondary coagulation standard that is traceable to current WHO International Standards.

B. Reagent Preparation/Handling:

Thaw each vial at 37 °C (+/- 1 °C) in a waterbath for 5 minutes. Invert gently before use.

C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60 °C or lower.

Once thawed, the plasma is stable in the original vial for 8 hours at 2-8 °C, except Protein S which is stable for 4 hours at 2-8 °C. Thawed material should be discarded after use and not be refrozen.

D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found nonreactive for Hepatitis B surface antigen (Hb_sAg), negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA) as well as for Hepatitis C (HCV). As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

Procedure

A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay Reagents

Quality Control Material (e.g. VisuCon-F Frozen Normal Control Plasma, VisuCon-F Frozen Abnormal Control Plasma)

Common clinical laboratory equipment and material (pipettes, etc.)

B. Assay Procedure

After thawing the Calibrator Plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established assay procedures for the quantitative determination of coagulation parameters.

Expected Values

Refer to the lot specific Certificate of Analysis for the assigned values for each parameter.

Quality Control

Each laboratory should establish and maintain its own quality control ranges for the particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, control plasmas or instrumentation to try and identify potential sources of error. If necessary, a new reference curve may need to be established.

Limitations of the Procedure

The results obtained with the VisuCal-F Frozen Calibrator Plasma will be subject to the limitations of the assay system including types of reagents, factor deficient substrates and instrumentation. If the lot of calibrator plasma or other critical reagent is changed, a new reference curve must be established.

Performance Characteristics

Values stated on the Certificate of Analysis are lot specific and vary from one lot to another.

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Bibliography/ Bibliographie / Bibliografía/ Bibliographie/ Bibliografía

1. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
2. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2012, European Committee for Standardization, 2012.

Limited Warranty: This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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AFFINITY BIOLOGICALS INC.

1395 Sandhill Drive
Ancaster ON CANADA L9G 4V5

Tel: (905) 304-9896
(800) 903-6020

Fax: (905) 304-9897
info@affinitybiologicals.com

EC

REP

Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands