# VisuCon™-L Coag Screen N

**Product #:** LCSN01

Store at 2-8°C.

 1348 Sandhill Drive. Ancaster, Ontario, Canada L9C 4V5  
 905-304-9896 • 800-903-6020 • fax 905-304-9897

 For *in vitro* diagnostic use.

## Intended Use

ENGLISH

The VisuCon-L Coag Screen N plasma is a normal control plasma intended for use in the quality control of quantitative coagulation assays, including Prothrombin Time (PT), Activated partial thromboplastin time (APTT) and fibrinogen, in the normal range. The VisuCon-L Coag Screen N plasma may be used with mechanical and photo-optical coagulation instruments in conjunction with appropriate commercial reagents.

## Summary and Principle

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-L Coag Screen N plasma is a pool of normal citrated human plasma collected from a minimum of 20 donors, buffered with 0.05M HEPES buffer, as well as stabilizing agents, dispensed and lyophilized. This control plasma is intended for use in the monitoring of the performance of PT, APTT and fibrinogen assays in the normal range.

## Reagents

### A. Description:

*Coag Screen N:*

Box containing 10 x 1mL vials of lyophilized plasma, Prod. # LCSN01

### B. Reagent Preparation/Handling:

Reconstitute each vial with 1 ml reagent grade water. Swirl gently and let stand for 15 minutes. Swirl gently once more before use.

### C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at 2-8°C. Once reconstituted, the plasma is stable for 24 hours on-board (15 - 22°C) or at 2-8°C in original vial. Gently mix vial of plasma if left for more than 2 hours without use. Reconstituted material should be discarded after use and not be refrozen.

### D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

## Procedure

### A. Materials Required (but not provided)

Coagulation Instrument or assay system

Assay Reagents

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Common clinical laboratory equipment and material (pipettes, etc.)

### B. Assay Procedure

After reconstituting the VisuCon-L Coag Screen N plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established coagulation assay procedures.

## Quality Control

Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to identify potential sources of error. QC materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

## Limitations of the Procedure

The results obtained with the VisuCon-L Coag Screen N plasma will be subject to the limitations of the assay system including types of reagents and instrumentation.

## Expected Values

The results obtained for the VisuCon-L Coag Screen N plasma will depend on various factors including assay method, types of reagents and instrumentation. Therefore, each laboratory should establish its own mean values and expected reference ranges for quality control purposes using their particular instrument-reagent system. The results presented in the supplied Certificate of Analysis were determined from data specific to the instrument-reagent system used for internal quality control testing. Lot specific results may vary when using different instrument and reagent combinations.

## Performance Characteristics

### PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuCon-L Coag Screen N plasma using the following instruments and reagents: STA® Neoplastine Cl+, HemosIL™ APTT-SP and STA® Fibrinogen-5 on two STA Compacts and Siemens Innovin and Actin FSL on two BCS systems. Each lot of VisuCon-L Coag Screen N was tested using the indicated instruments and reagents in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2<sup>2</sup>. The coefficients of variation (% CV) obtained in these precision studies are presented in the tables below.

### STA Compact, Neoplastine Cl+, APTT-SP and Fibrinogen 5

	Within Run	Between Day	Between Run	Within Device
PT	0.8 - 1.5%	0.5 - 1.6%	0.3 - 1.0%	1.7 - 2.0%
APTT	1.0 - 1.8%	0.3 - 1.6%	0 - 1.9%	1.8 - 2.3%
Fibrinogen	2.0 - 3.0%	0 - 4.0%	1.0 - 1.9%	3.0 - 5.1%

### BCS, Innovin and Actin FSL

	Within Run	Between Day	Between Run	Within Device
PT	1.0 - 1.5%	0.6 - 1.1%	0.8 - 1.3%	1.5 - 2.1%
APTT	0.6 - 1.1%	0 - 0.6%	0.4 - 1.0%	0.9 - 1.5%

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---

1. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.
2. "Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition". EP5-A2, CLSI, Vol. 24, No. 25, 2004.

**SYMBOL LEGEND / LEGENDE DES SYMBOLES / LEGENDA DEI  
SIMBOLI / SYMBOL LEGENDE / LEYENDAS DE SÍMBOLOS /  
VERKLARING VAN DE SYMBOLEN**<sup>1</sup>

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

Temperature limitation / Limite de température / Limite di temperatura / Temperaturgrenzwert / Limite de temperature / Temperatuurbeperring

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

**Limited Warranty:**

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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