****REPRESENTATIVE DATASHEET****

VisuCon™-F Frozen Unassayed Normal Control Plasma

Product #: UFNCP0125
UFNCP0181
UFNCP0481

Store at -60°C or lower.

For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuCon-F Frozen Unassayed Normal Control plasma is intended for use in the quality control of coagulation assays in the normal range.

Summary and Principle

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-F Frozen Unassayed Normal Control plasma is a pool of normal citrated human plasma collected from a minimum of 20 donors, buffered with 0.02M HEPES buffer, dispensed and rapidly frozen. This control plasma is intended for use in the monitoring of the performance of PT, APTT and fibrinogen assays in the normal range.

Reagents

A. Description:

Unassayed Normal Control Plasma:

Box containing 25 x 1mL vials of frozen plasma, Prod. # UFNCP0125
OR
Box containing 81 x 1mL vials of frozen plasma, Prod. # UFNCP0181
OR
Box containing 81 x 4mL vials of frozen plasma, Prod. # UFNCP0481

B. Reagent Preparation/Handling:

Thaw each 1 mL vial for 5 minutes and each 4 mL vial for 7 minutes in a 37°C (+/- 1°C) waterbath. Invert gently before use.

C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower.

Once thawed, the plasma is stable for 8 hours on-board or at 2-8°C in original vial. Thawed material should be discarded after use and not be refrozen.

D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA) as well as for Hepatitis C (HCV). As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

Procedure

A. Materials Required (but not provided)

Waterbath
Coagulation Instrument or assay system
Assay Reagents
Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)
Common clinical laboratory equipment and material (pipettes, etc.)

B. Assay Procedure

After thawing the VisuCon-F Frozen Unassayed Normal Control Plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established coagulation assay procedures.

Expected Values

The results obtained for the unassayed normal control plasma will depend on various factors including types of reagents, factor deficient substrates and instrumentation therefore each laboratory should establish its own mean values and expected reference ranges for quality control purposes in every coagulation assay using their particular instrument-reagent system.

Quality Control

Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Limitations of the Procedure

The results obtained with the VisuCon-F Frozen Unassayed Normal Control plasma will be subject to the limitations of the assay system including types of reagents, factor deficient substrates and instrumentation.

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Bibliography/Bibliographie/

Bibliografia/Bibliographie/ Bibliografia

1. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
2. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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