**\*\*REPRESENTATIVE DATASHEET\*\***

**VisuCon™-F Frozen Abnormal Control Plasma**

**Product #:** FRABP0125  
FRABP0105

Store at -60°C or lower.

For *In Vitro* Diagnostic Use.

**Intended Use**

ENGLISH

The VisuCon-F Frozen Abnormal Control plasma is intended for use in the quality control of coagulation assays in the borderline pathological range.

**Summary and Principle**

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-F Frozen Abnormal Control plasma is a pool of a minimum of 20 normal citrated human plasma donors, diluted to defined concentrations, buffered with 0.02M HEPES buffer and rapidly frozen. This control plasma may be used to monitor the performance of coagulation assays for the following parameters: fibrinogen (Clauss Method), coagulation factors II, V, VII, VIII, IX, X, XI, XII, antithrombin, Protein C, and Protein S.

**Reagents**

**A. Description:**

*Abnormal Control Plasma:*

Box containing 25 x 1mL vials of frozen plasma, Prod. # FRABP0125

OR

Box containing 5 x 1mL vials of frozen plasma, Prod. # FRABP0105

Each box of VisuCon-F Frozen Abnormal Control Plasma contains a lot-specific Certificate of Analysis with the assigned reference ranges for each parameter. The reported reference ranges were determined against, where available, the SSC/ISTH secondary coagulation standard that is traceable to current WHO International Standards.

**B. Reagent Preparation/Handling:**

Thaw each vial at 37°C (+/- 1°C) in a waterbath for 5 minutes. Invert gently before use.

**C. Storage and Stability:**

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower.

Once thawed, the plasma is stable for 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and not be refrozen.

**D. Precautions and Warnings:**

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA) as well as for Hepatitis C (HCV). As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

**Procedure**

**A. Materials Required (but not provided)**

Waterbath

Coagulation Instrument or assay system

Assay Reagents

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Common clinical laboratory equipment and material (pipettes, etc.)

**B. Assay Procedure**

After thawing the Control Plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established assay procedures for the quantitative determination of coagulation parameters.

**Expected Values**

Refer to the lot specific Certificate of Analysis for the assigned range for each parameter. The results obtained for the control plasma will depend on various factors including types of reagents, factor deficient substrates and instrumentation therefore it is recommended that each laboratory establish its own mean values and expected reference ranges for quality control purposes using their particular instrument-reagent system.

**Quality Control**

Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

**Limitations of the Procedure**

The results obtained with the VisuCon-F Frozen Abnormal Control plasma will be subject to the limitations of the assay system including types of reagents, factor deficient substrates and instrumentation.

**Performance Characteristics**

Values stated on the Certificate of Analysis are lot specific and vary from one lot to another.

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## Bibliography / Bibliographie / Bibliografia / Bibliographie / Bibliografía

1. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
2. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

## SYMBOL LEGEND/ LEGENDE DES SYMBOLES/ LEGENDA DEI SIMBOLI/ SYMBOL LEGENDE/ LEYENDAS DE SÍMBOLOS 2

 For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro

 Batch code / Code du lot / Codice lotto / Chargencode / Código de lote

 Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado

 Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad

 Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo

 Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenwert / Límite superior de temperatura

 Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante

 Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso

 Biological Risks / Risques biologiques / Rischio Biológico / Biologische Risiken / Riesgos biológicos

**Limited Warranty:** This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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