**\*\*REPRESENTATIVE DATASHEET\*\***
**VisuCon™-F Low Fibrinogen Control Plasma**
**Product #:** LFGCP0125  
 LFGCP0105  
 LFGCP0181

Store at -60°C or lower.

 For *in vitro* diagnostic use.

**Intended Use**

ENGLISH

The VisuCon-F Low Fibrinogen Control Plasma is an assayed control plasma prepared from de-fibrinated human plasma intended for use in the quality control of quantitative fibrinogen assays in the low abnormal range. The VisuCon-F Low Fibrinogen Control Plasma may be used with mechanical instruments in conjunction with appropriate commercial reagents for determining fibrinogen levels in plasma by the clotting method of Clauss. This plasma is intended "For *In Vitro* Diagnostic Use".

The intended users of the VisuCon-F Low Fibrinogen Control Plasma are trained laboratory personnel working in clinical laboratories.

**Summary and Principle**

Fibrinogen is an abundant plasma protein (5-10 µM) produced in the liver with a molecular weight of 340 kDa. Proteolysis of fibrinogen by thrombin results in the release of fibrinopeptide A followed by fibrinopeptide B and the formation of fibrin monomers that polymerize to form insoluble fibrin fibrils. The polymerised fibrin is subsequently stabilized by activated Factor XIII. Fibrin formation and its subsequent degradation are carefully regulated processes; therefore, changes in either may lead to hemorrhage or thrombosis. As fibrinogen is an acute phase reactant stresses including trauma, pregnancy and inflammation will result in an increase in its concentration. Elevated levels can also occur with contraceptive use and in cardiovascular disease. Congenital disorders of fibrinogen are rare. Acquired disorders are fairly common in liver disease and decreased levels may also be observed in DIC.<sup>1-4</sup>

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-F Low Fibrinogen Control plasma is a pool of de-fibrinated citrated human plasma buffered with 0.02M HEPES buffer, dispensed and rapidly frozen. This control plasma is used to monitor the performance of fibrinogen assays in the low abnormal range.

**Reagents**
**A. Description:**
*Low Fibrinogen Control Plasma:*

 Box containing 25 x 1 mL vials of frozen plasma, Prod. # LFGCP0125  
 OR

 Box containing 5 x 1 mL vials of frozen plasma, Prod. # LFGCP0105  
 OR

Box containing 81 x 1 mL vials of frozen plasma, Prod. # LFGCP0181

**B. Reagent Preparation/Handling:**

Thaw each 1 ml vial at 37°C (+/- 1°C) in a waterbath for 5 minutes. Invert gently before use.

**C. Storage and Stability:**

The VisuCon-F Low Fibrinogen Control Plasma is shipped frozen on dry ice and should be stored immediately at -60°C or lower upon receipt. If plasma is not received in a frozen state, it should be discarded and not used for quality control testing.

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, the plasma is stable for 8 hours on-board (19-22°C) or for 72 hours at 2-8°C in the original vial. Thawed material should be discarded after use and not be refrozen.

**D. Precautions and Warnings:**

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found

non-reactive for Hepatitis B surface antigen (HbsAg), negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA) as well as for Hepatitis C (HCV). As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

**Procedure**
**A. Materials Required (but not provided)**

 Waterbath  
 Diagnostica Stago STA Compact Analyzer  
 STA Fibrinogen 5 Assay reagents  
 Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)  
 Common clinical laboratory equipment and material (pipettes, etc.)

**B. Assay Procedure**

After thawing, the VisuCon-F Low Fibrinogen Control Plasma can be used on the Diagnostica Stago STA Compact analyzer using Diagnostica Stago STA Fibrinogen-5 reagent in accordance with established assay procedures for the quantitative determination of fibrinogen.

**Expected Values**

Each box of VisuCon-F Low Fibrinogen Control Plasma contains a lot-specific Certificate of Analysis with an assigned range for fibrinogen. The reported range was determined against the SSC/ISTH secondary coagulation standard that is traceable to current WHO International Standards. The results obtained for the control plasma will depend on various factors including types of reagents and instrumentation, therefore, it is recommended that each laboratory establish its own mean values and reference ranges for quality control purposes using their particular instrument-reagent system.

**Quality Control**

Each laboratory should establish and maintain its own quality control ranges for the particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, calibrator plasma or instrumentation to identify potential sources of error.

**Limitations of the Procedure**

The results obtained with the VisuCon-F Low Fibrinogen Control plasma will be subject to the limitations of the assay system including calibrator plasma, types of reagents and instrumentation.

**Performance Characteristics**
**PRECISION:**

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuCon-F Low Fibrinogen Control Plasma using Diagnostica Stago Fib-5 Reagent on 2 STA Compacts. Each lot of plasma was tested in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2<sup>5</sup>. The coefficients of variation (% CV) obtained in this precision study are presented in the table below.

Within Run	Between Day	Between Run	Within Device
3.15%	0.54%	0%	8.44%

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1. Hantgan RR, Francis CW, Marder VJ; Fibrinogen Structure and Physiology; in Hemostasis and Thrombosis, 3<sup>rd</sup> Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp 277-300, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
2. Shafer JA, Higgins DL; Human Fibrinogen; CRC Critical Reviews in Clinical Laboratory Sciences 26, pp 1-41, 1988.
3. Binnie CG, Lord ST; The Fibrinogen Sequences that Interact with Thrombin; Blood 81, pp 3186-3192, 1993.
4. McDonagh J, Carrell N, Lee MH; Dysfibrinogenemia and Other Disorders of Fibrinogen Structure and Function; in Hemostasis and Thrombosis, 3<sup>rd</sup> Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp 314-334, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
5. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
6. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

---

**Limited Warranty:** This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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