****REPRESENTATIVE DATASHEET****

VisuCon™-F Frozen Normal Control Plasma

Product #: FRNCP0125
FRNCP0105

Store at -60°C or lower.

For *In Vitro* Diagnostic Use.

Intended Use

ENGLISH

The VisuCon-F Frozen Normal Control plasma is intended for use in the quality control of coagulation assays in the normal range.

Summary and Principle

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-F Frozen Normal Control plasma is a pool of normal citrated human plasma collected from a minimum of 20 donors, buffered with 0.02M HEPES buffer, dispensed and rapidly frozen. This control plasma may be used to monitor the performance of coagulation assays for the following parameters: fibrinogen (Claus Method), coagulation factors II, V, VII, VIII, IX, X, XI, XII, antithrombin, Protein C, and Protein S.

Reagents

A. Description:

Normal Control Plasma:

Box containing 25 x 1mL vials of frozen plasma, Prod. # FRNCP0125

OR

Box containing 5 x 1mL vials of frozen plasma, Prod. # FRNCP0105

Each box of VisuCon-F Frozen Normal Control Plasma contains a lot-specific Certificate of Analysis with the assigned reference ranges for each parameter. The reported reference ranges were determined against, where available, the SSC/ISTH secondary coagulation standard that is traceable to current WHO International Standards.

B. Reagent Preparation/Handling:

Thaw each vial at 37°C (+/- 1°C) in a waterbath for 5 minutes. Invert gently before use.

C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower.

Once thawed, the plasma is stable for 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and not be refrozen.

D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA) as well as for Hepatitis C (HCV). As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

Procedure

A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay Reagents

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Common clinical laboratory equipment and material (pipettes, etc.)

B. Assay Procedure

After thawing the Control Plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established assay procedures for the quantitative determination of coagulation parameters.

Expected Values

Refer to the lot specific Certificate of Analysis for the assigned range for each parameter. The results obtained for the control plasma will depend on various factors including types of reagents, factor deficient substrates and instrumentation therefore it is recommended

that each laboratory establish its own mean values and expected reference ranges for quality control purposes using their particular instrument-reagent system.

Quality Control

Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Limitations of the Procedure

The results obtained with the VisuCon-F Frozen Normal Control plasma will be subject to the limitations of the assay system including types of reagents, factor deficient substrates and instrumentation.

Performance Characteristics

Values stated on the Certificate of Analysis are lot specific and vary from one lot to another.

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Bibliography/ Bibliographie / Bibliografia / Bibliographie / Bibliografía

1. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
2. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

SYMBOL LEGEND/ LEGENDE DES SYMBOLES/ LEGENDA DEI SIMBOLI / SYMBOL LEGENDE /LEYENDAS DE SÍMBOLOS 2

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo

Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenwert / Limite superior de temperatura

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos

fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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Ermgo Europe

Prinsessegracht 20
2514 AP The Hague
The Netherlands

Limited Warranty: This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or