VisuCon™-L Abnormal Control Plasma

Product #: LYABP01

1348 Sandhill Drive. Ancaster, Ontario, Canada L9G 4V5
905-304-9896 • 800-903-6020 • fax 905-304-9897

Store at 2-8°C.

For *In Vitro* Diagnostic Use.

Intended Use

ENGLISH

The VisuCon-L Abnormal Control plasma is an assayed control plasma intended for use in the quality control of coagulation assays in the borderline pathological range. The VisuCon-L Abnormal Control Plasma may be used with mechanical or photo-optical instruments in conjunction with appropriate commercial reagents for measuring fibrinogen (Clauss Method), coagulation factors II, V, VII, VIII, IX, X, XI, XII, antithrombin activity, Protein C activity, and Protein S activity. This plasma is intended "For *In Vitro* Diagnostic Use".

The intended users of the VisuCon-L Abnormal Control Plasma are trained laboratory personnel working in clinical laboratories.

Summary and Principle

Good laboratory practice requires that assays be calibrated over the relevant range and that controls be performed regularly to confirm assay calibration. The VisuCon-L Abnormal Control plasma is a pool of a minimum of 20 normal citrated human plasma donors, diluted to defined concentrations, buffered with 0.05M HEPES buffer, as well as stabilizing agents, dispensed and lyophilized. This control plasma may be used to monitor the performance of quantitative coagulation assays for the following parameters: fibrinogen (Clauss Method), coagulation factors II, V, VII, VIII, IX, X, XI, XII, antithrombin activity, Protein C activity, and Protein S activity.

Reagents

A. Description:

Abnormal Control Plasma:

Box containing 10 x 1mL vials of lyophilized plasma, Prod. # LYABP01

Each box of VisuCon-L Abnormal Control Plasma contains a lot-specific Certificate of Analysis with the assigned reference ranges for each parameter. The reported reference ranges were determined against, where available, the SSC/ISTH secondary coagulation standard, that is traceable to current WHO International Standards or an approved diagnostic calibrator.

B. Reagent Preparation/Handling:

Reconstitute each vial with 1ml reagent grade water. Swirl gently and let stand for 15 minutes. Swirl gently once more before use.

C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at 2-8°C. Once reconstituted, the plasma is stable for 8 hours at 2-8°C or on-board (15-22°C) in original vial, except Protein S, which is stable for 4 hours at both temperatures. Gently mix vial

of plasma if left for more than 2 hours without use. Reconstituted material should be discarded after use and not be refrozen.

D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

Procedure

A. Materials Required (but not provided)

Reagent Grade Water

Coagulation Instrument or assay system

Assay Reagents

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Common clinical laboratory equipment and material (pipettes, etc.)

B. Assay Procedure

After reconstituting the VisuCon-L Abnormal Control Plasma, use on coagulation instrument or assay system with corresponding reagents in accordance with established assay procedures for the quantitative determination of coagulation parameters.

Expected Values

The results obtained for the VisuCon-L Abnormal Control plasma will depend on various factors including types of reagents, factor deficient substrates and instrumentation. Therefore, each laboratory should establish its own mean values and expected reference ranges for quality control purposes using their particular instrument-reagent system. The results presented in the supplied Certificate of Analysis were determined from data specific to the instrument-reagent system used for internal quality control testing. Lot specific results may vary when using different instrument and reagent combinations.

Quality Control

Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error. QC materials should be used in accordance with local, state and/or federal regulations or accreditation requirements.

Limitations of the Procedure

The results obtained with the VisuCon-L Abnormal Control plasma will be subject to the limitations of the assay system including types of reagents, factor deficient substrates and instrumentation.

Performance Characteristics

Values stated on the Certificate of Analysis are lot specific and vary from one lot to another.

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1. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

Symbol Legend/ Legende des Symboles/ Legenda dei Simboli / Symbol Legende / Leyendas de Simbolos / Verklaring van de Symbolen¹

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

Temperature limitation / Limite de température / Limite di temperatura / Temperaturgrenzwert / Limite de temperature / Temperatuurbepanking

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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