****REPRESENTATIVE DATASHEET****
VisuDep™-F Frozen Factor IX Deficient Plasma
Product #: FRDP090125
FRDP090105

Store at -60°C or lower.

 For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuDep-F Frozen Factor IX Deficient Plasma is intended for use as a substrate plasma in the quantitative determination of Factor IX activity in citrated plasmas using a clot-based assay.

Summary and Principle

 Factor IX (Christmas Factor) is a vitamin K-dependent glycoprotein that circulates in plasma at a concentration of 5 µg/ml.¹ Factor IX contains two EGF-like domains and an amino-terminal domain containing 12 γ-carboxy-glutamic acid (Gla) residues. These Gla residues allow Factor IX to bind divalent metal ions and participate in calcium-dependent binding interactions.¹ Factor IX can be activated by limited proteolysis in the presence of calcium by activated factor XI and/or through the extrinsic coagulation pathway by a complex of VII^a/tissue factor/phospholipid and activated Factor X.^{1,2} Activated Factor IX, along with Factor VIIIa, calcium ions and a phospholipid membrane, converts Factor X into Xa eventually leading to the formation of a fibrin clot.³

 The biological importance of Factor IX is demonstrated in Haemophilia B (Christmas disease), an X-linked congenital bleeding disease resulting from a quantitative (low activity and low antigen) or qualitative (low activity and normal antigen) defect in Factor IX function.⁴ The congenital deficiency of Factor IX may be classified as severe (<1% Factor IX activity), moderate (between 1 and 5% Factor IX activity) or mild (between 5 and 40% Factor IX activity).⁵

 Factor IX activity in patient plasma is determined *in vitro* through the use of a modified activated partial thromboplastin time test (APTT). Patient plasma is diluted, mixed with plasma deficient in Factor IX and an APTT test performed. The result is interpolated using a reference curve generated with dilutions of calibrator plasma mixed with the Factor IX deficient plasma. Correction of the clotting time of the deficient plasma is proportional to the activity level of Factor IX in the patient plasma.

Reagents
A. Description:

Frozen Factor IX immune-depleted plasma with residual Factor IX activity of <1%. Remaining coagulation factors have activities of ≥50% and fibrinogen is ≥1.5 g/L.

Box containing 25 x 1 ml vials of frozen Factor IX deficient plasma,

Product # FRDP090125

OR

Box containing 5 x 1 ml vials of frozen Factor IX deficient plasma,

Product # FRDP090105
B. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that

products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

C. Reagent Preparation/Handling:

Thaw each 1 mL vial for 5 minutes in a 37°C ± 1°C waterbath. Invert gently before use.

D. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma is stable for 4 hours on-board (15-19°C) or 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and should not be refrozen.

Specimen Collection and Preparation

 Blood is collected into 3.2% sodium Citrate anticoagulant tubes at a ratio of 9 volumes blood to 1 volume anticoagulant and gently mixed by inversion. Refer to CLSI Guideline H21-A5 "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" for further information on specimen collection, handling and storage.⁶
Procedure
A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay reagents e.g. APTT reagent, calcium chloride, Owren's buffer

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Control Plasmas (e.g. VisuCon-F Frozen Normal and Abnormal Control Plasmas)

Common clinical laboratory equipment and material (pipettes, plastic tubes etc.)

B. Assay Procedure (Manual Method)
Preparation of Calibration Curve Dilutions

i) Prepare calibrator plasma, APTT reagent and assay buffer according to manufacturer's instructions.

ii) Prepare dilutions of the calibrator plasma in plastic tubes using assay buffer to generate a 5-point calibration curve corresponding to 100% to 6.25%, as illustrated in the following table:

Dilution	% Factor IX *	Calibrator Plasma Volume	Buffer Volume
1/10	100%	100 µl	900 µl
1/20	50%	500 µl of 1/10 dilution	500 µl
1/40	25%	500 µl of 1/20 dilution	500 µl
1/80	12.5%	500 µl of 1/40 dilution	500 µl
1/160	6.25%	500 µl of 1/80 dilution	500 µl

*Refer to manufacturer's Certificate of Analysis for the lot-specific Factor IX activity level of the calibrator plasma.

Preparation of Test Plasma Dilutions

 Dilute all test plasma samples 1/10 and 1/20 in plastic tubes using assay buffer. Additional dilutions may be required if the reported results are outside the readable range of the calibration curve or if the test results are not parallel with the calibration curve.⁷
Test Method (Manual)

All calibration plasma and test plasma dilutions should be tested immediately after preparation. Add the following to a test tube or cuvette warmed to 37°C ± 1°C:

- 100 µl of VisuDep-F Factor IX deficient plasma
- 100 µl of calibrator plasma dilution or diluted test sample

Mix and incubate for 1 minute at 37°C ± 1°C then add:

- 100 µl of APTT reagent (pre-warmed to 37°C ± 1°C)

Mix and incubate for 3-5 minutes at 37°C ± 1°C then add:

- 100 µl of calcium chloride (pre-warmed to 37 °C ± 1 °C)

Start a timing device the instant the calcium chloride is added to the plasma mixture and record the time to clot formation.

Results

Using log-log graph paper, plot clot times for the calibrator plasma dilutions on the y-axis against the corresponding percent Factor IX activity on the x-axis. Draw the best-fit calibration line. Interpolate the Factor IX levels of the test samples from this calibration curve and multiply by the appropriate dilution factor. Under the conditions described here, a test sample diluted 1/10 will have a dilution factor of 1 and a dilution of 1/20 will have a dilution factor of 2. Additionally, if the Factor IX value for the calibrator plasma is not 100%, the test sample result will need to be multiplied by the appropriate correction factor. For example, if the Factor IX level of the calibrator plasma is 105%, multiply the test sample results by 1.05 to obtain the corrected result.

Automatic Test Method

Factor activity assays may be performed using semi-automated or automated coagulation instruments with built in programs for generating the calibration curve and calculating test plasma results. Refer to the coagulation instrument's Operator's Manual for specific assay procedure instructions.

Quality Control

In order to ensure accuracy of the results, normal and abnormal control plasmas such as the VisuCon-F Frozen Normal and Abnormal Controls, should be analyzed regularly. Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Results

Results for test samples and control plasmas can be calculated as described above or using the software program of the coagulation instrument. Ensure values obtained for the control plasmas are within established ranges. If control results are out of range, the Factor IX activity values obtained for test samples should be considered suspect.

Limitations of the Procedure

Incorrect sample handling such as insufficient filling of collection tubes, contamination of collection tubes or incorrect type of anticoagulant can lead to erroneous sample results. Use only plasma samples that have been collected into sodium citrate anticoagulant.⁶ Samples that have clotted or are hemolyzed, icteric or lipemic should be avoided. Do not use glass tubes to prepare plasma dilutions as glass may cause activation of the coagulation cascade.⁷⁻⁸ The presence of lupus anticoagulant in test plasmas can affect the reported factor activity.

Expected Values

The normal range for Factor IX as reported in the literature is typically 55-163%⁹, however, each laboratory should determine a normal range independently. Results measured in 100 healthy individuals indicate a Factor IX normal reference range of 75 - 159%.

Performance Characteristics

PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuDep-F Frozen Factor IX Deficient Plasma using STA® PTT-Automate on an STA Compact. Normal and abnormal control plasmas were tested in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2¹⁰. The coefficients of variation (% CV) obtained in these precision studies are presented in the table below.

	Within Run	Between Day	Between Run	Within Device
Normal	6.0 - 7.4%	5.1 - 6.4%	4.1 - 6.8%	8.8 - 11.9%
Abnormal	4.0 - 5.1%	5.2 - 5.9%	7.0 - 7.8%	9.6 - 10.7%

METHOD COMPARISON:

In a method comparison study, 120 test samples containing Factor IX ranging from <1% to normal levels were tested using the VisuDep-F Frozen Factor IX Deficient Plasma and Precision Biologic's Factor IX Deficient Plasma as the reference method. The correlation co-efficient (r) reported in this study was 0.991 (R² = 0.982, y = 0.956x + 2.589).

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
1. Limentani SA, Furie BC, Furie B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 94-108, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
2. Lawson JH, Mann KG; Cooperative Activation of Human Factor IX by the Human Extrinsic Pathway of Coagulation; JBC 266 pp11317-11327, 1991.
3. van Hylckama Vlieg A, van der Linden IK, Bertina RM, Rosendaal FR. High levels of Factor IX increase the risk of venous thrombosis, *Blood*, 2000, 95:12, pp. 3678-3682.
4. Brettler DB, Levine PH, "Clinical Manifestations and Therapy of Inherited Coagulation Factor Deficiencies" in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 169-183, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
5. White GC, Rosendaal F, Aledort LM, Lusher JM, Rothschild C, Ingerslev J. Definitions in Hemophilia, Recommendation of the Scientific Subcommittee on Factor VIII and Factor IX of the Scientific and Standardization Committee of the International Society on Thrombosis and Haemostasis, *Thrombosis and Haemostasis*, 2001, 85, p. 560.
6. Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General performance of Coagulation Assays; Approved Guideline—Fifth Edition. CLSI Document H21-A5, Vol. 28, No. 5. 2008.
7. "Determination of Coagulation Factor Activities Using the One-Stage Clotting Assay", H48, 2nd edition, CLSI, March 2016.
8. "One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test", Approved Guideline, Second Edition, H47-A2, CLSI, Vol. 28, No. 20, 2008.
9. Andrew, M, Schmidt, B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 991-992, J.B. Lippincott Co., Philadelphia, 1994.
10. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
11. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

**SYMBOL LEGEND/ LEGENDE DES SYMBOLES/ LEGENDA DEI
SIMBOLI/ SYMBOL LEGENDE/ LEYENDAS DE SÍMBOLOS /
VERKLARING VAN DE SYMBOLEN 11**

 For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

 Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

 Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

 Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

 Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

 Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenwert / Limite superior de temperatura / Bovenste temperatuurbegrenzing

 Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

 Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

 Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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