****REPRESENTATIVE DATASHEET****
VisuDep™-F Frozen Factor VIII Deficient Plasma
Product #: FRDP080125
FRDP080105

Store at -60°C or lower.

 1348 Sandhill Drive, Ancaster, Ontario, Canada L9C 4V5
905-304-9896 • 800-903-6020 • fax 905-304-9897

 For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuDep-F Frozen Factor VIII Deficient Plasma is intended for use as a substrate plasma in the quantitative determination of Factor VIII activity in citrated plasmas using a clot-based assay.

Summary and Principle

Factor VIII is a large glycoprotein with a molecular weight of 320000 daltons that circulates in plasma at a concentration of approximately 200 ng/mL.^{1,2} FVIII is stabilized by association with von Willebrand Factor to form a FVIII-vWF complex required for the normal survival of FVIII *in vivo* ($t_{1/2}$ of 8-12 hours).³ FVIII is a pro-cofactor that is activated through limited proteolysis by thrombin. In this process FVIIIa dissociates from vWF to combine with activated Factor IX, calcium and a phospholipid surface where it is an essential cofactor in the assembly of the Factor X activator complex.^{4,5} Once dissociated from vWF, FVIIIa is susceptible to inactivation by activated Protein C and by non-enzymatic decay.⁴

The biological importance of Factor VIII is demonstrated in Hemophilia A, a congenital bleeding disorder occurring primarily in males that results from an X-chromosome-linked deficiency of FVIII.⁶ The severity of the deficiency generally correlates with the severity of the disease. Individuals with <1% Factor VIII activity are classified as severe patients, those with between 1 and 5% Factor VIII activity are classified as moderate and those with between 5 and 40% Factor VIII activity are classified as mild hemophiliacs.⁷

Factor VIII activity in patient plasma is determined *in vitro* through the use of a modified activated partial thromboplastin time test (APTT). Patient plasma is diluted, mixed with plasma deficient in Factor VIII and an APTT test performed. The result is interpolated using a reference curve generated with dilutions of calibrator plasma mixed with the Factor VIII deficient plasma. Correction of the clotting time of the deficient plasma is proportional to the activity level of Factor VIII in the patient plasma.

Reagents
A. Description:

 Frozen Factor VIII immune-depleted plasma with residual Factor VIII activity of <1%. Remaining coagulation factors have activities of $\geq 50\%$ and fibrinogen is ≥ 1.5 g/L.

Box containing 25 x 1 ml vials of frozen Factor VIII deficient plasma,

Product # FRDP080125

OR

Box containing 5 x 1 ml vials of frozen Factor VIII deficient plasma,

Product # FRDP080105
B. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

C. Reagent Preparation/Handling:

 Thaw each 1 mL vial for 5 minutes in a 37°C \pm 1°C waterbath. Invert gently before use.

D. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma is stable for 4 hours on-board (15-19°C) or 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and should not be refrozen.

Specimen Collection and Preparation

Blood is collected into 3.2% sodium Citrate anticoagulant tubes at a ratio of 9 volumes blood to 1 volume anticoagulant and gently mixed by inversion. Refer to CLSI Guideline H21-A5 "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" for further information on specimen collection, handling and storage.⁸

Procedure
A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay reagents e.g. APTT reagent, calcium chloride, Owren's buffer

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Control Plasmas (e.g. VisuCon-F Frozen Normal and Abnormal Control Plasmas)

Common clinical laboratory equipment and material (pipettes, plastic tubes etc.)

B. Assay Procedure (Manual Method)
Preparation of Calibration Curve Dilutions

i) Prepare calibrator plasma, APTT reagent and assay buffer according to manufacturer's instructions.

ii) Prepare dilutions of the calibrator plasma in plastic tubes using assay buffer to generate a 5-point calibration curve corresponding to 100% to 6.25%, as illustrated in the following table:

Dilution	% Factor VIII *	Calibrator Plasma Volume	Buffer Volume
1/10	100%	100 μ l	900 μ l
1/20	50%	500 μ l of 1/10 dilution	500 μ l
1/40	25%	500 μ l of 1/20 dilution	500 μ l
1/80	12.5%	500 μ l of 1/40 dilution	500 μ l
1/160	6.25%	500 μ l of 1/80 dilution	500 μ l

*Refer to manufacturer's Certificate of Analysis for the lot-specific Factor VIII activity level of the calibrator plasma.

Preparation of Test Plasma Dilutions

Dilute all test plasma samples 1/10 and 1/20 in plastic tubes using assay buffer. Additional dilutions may be required if the reported results are outside the readable range of the calibration curve or if the test results are not parallel with the calibration curve.⁹

Test Method (Manual)

All calibration plasma and test plasma dilutions should be tested immediately after preparation. Add the following to a test tube or cuvette warmed to 37°C \pm 1°C:

- 100 μ l of VisuDep-F Factor VIII deficient plasma
- 100 μ l of calibrator plasma dilution or diluted test sample

 Mix and incubate for 1 minute at 37°C \pm 1°C then add:

- 100 μ l of APTT reagent (pre-warmed to 37°C \pm 1°C)

 Mix and incubate for 3-5 minutes at 37°C \pm 1°C then add:

- 100 μ l of calcium chloride (pre-warmed to 37°C \pm 1°C)

Start a timing device the instant the calcium chloride is added to the plasma mixture and record the time to clot formation.

Results

Using log-log graph paper, plot clot times for the calibrator plasma dilutions on the y-axis against the corresponding percent Factor VIII activity on the x-axis. Draw the best-fit calibration line. Interpolate the Factor VIII levels of the test samples from this calibration curve and multiply by the appropriate dilution factor. Under the conditions described here, a test sample diluted 1/10 will have a dilution factor of 1 and a dilution of 1/20 will have a dilution factor of 2. Additionally, if the Factor VIII value for the calibrator plasma is not 100%, the test sample result will need to be multiplied by the appropriate correction factor. For example, if the Factor VIII level of the calibrator plasma is 105%, multiply the test sample results by 1.05 to obtain the corrected result.

Automatic Test Method

Factor activity assays may be performed using semi-automated or automated coagulation instruments with built in programs for generating the calibration curve and calculating test plasma results. Refer to the coagulation instrument's Operator's Manual for specific assay procedure instructions.

Quality Control

In order to ensure accuracy of the results, normal and abnormal control plasmas such as the VisuCon-F Frozen Normal and Abnormal Controls, should be analyzed regularly. Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Results

Results for test samples and control plasmas can be calculated as described above or using the software program of the coagulation instrument. Ensure values obtained for the control plasmas are within established ranges. If control results are out of range, the Factor VIII activity values obtained for test samples should be considered suspect.

Limitations of the Procedure

Incorrect sample handling such as insufficient filling of collection tubes, contamination of collection tubes or incorrect type of anticoagulant can lead to erroneous sample results. Use only plasma samples that have been collected into sodium citrate anticoagulant.⁸ Samples that have clotted or are hemolyzed, icteric or lipemic should be avoided. Do not use glass tubes to prepare plasma dilutions as glass may cause activation of the coagulation cascade.⁹⁻¹⁰ The presence of lupus anticoagulant in test plasmas can affect the reported factor activity.

Expected Values

The normal range for Factor VIII as reported in the literature is typically 50-149%¹¹, however, each laboratory should determine a normal range independently. Results measured in 94 healthy individuals indicate a Factor VIII normal reference range of 50-157%.

Performance Characteristics

PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuDep-F Frozen Factor VIII Deficient Plasma using STA® CK-Prest on an STA Compact. Normal and abnormal control plasmas were tested in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2¹². The coefficients of variation (% CV) obtained in these precision studies are presented in the table below.

	Within Run	Between Day	Between Run	Within Device
Normal	5.2 – 6.1%	0 – 5.2%	5.5 – 10.5%	9.2 – 12.2%
Abnormal	3.7 – 6.7%	0 – 4.9%	4.8 – 10.2%	7.8 – 12.2%

METHOD COMPARISON:

In a method comparison study, 150 test samples containing Factor VIII ranging from <1% to normal levels were tested using the VisuDep-F Frozen Factor VIII Deficient Plasma and Precision Biologic's Factor VIII Deficient Plasma as the reference method. The correlation co-efficient (r) reported in this study was 0.995 ($R^2 = 0.9893$, $y = 1.0223x + 0.7706$).

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1. Pittman, D.D., Kaufman, R.J. Structure-Function Relationships of Factor VIII Elucidated through Recombinant DNA Technology. *Thromb. Haemostas.*, 61(2), pp. 161-165, 1989.
2. Hoyer LW, Wyshock EG, Colman RW, "Coagulation Cofactors: Factor V and VIII" in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 109-133, J.B. Lippincott Co., Philadelphia, 1994.
3. Brettler, D.B., Levine, P.H., "Clinical Manifestations and Therapy of Inherited Coagulation Factor Deficiencies" in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 169-183, J. B. Lippincott Co., Philadelphia, 1994.
4. Lenting, P.J., van Mourik, J.A., Mertens, K. The Life Cycle of Coagulation Factor VIII in View of its Structure and Function. *Blood*, 92(11), pp. 3983-3996, 1998.
5. Furie B, Limentani SA, Rosenfield CG. A Practical Guide to the Evaluation and Treatment of Hemophilia, *Blood*, 1994, 84(1), pp. 3-9.
6. Hedner, U., Ginsburg, D., Lusher, J.M., High, K. A. Congenital Hemorrhagic Disorders: New Insights into the Pathophysiology and Treatment of Hemophilia. *Hematology*, pp. 241- 290, 2000.
7. White GC, Rosendaal F, Aledort LM, Lusher JM, Rothschild C, Ingerslev J. Definitions in Hemophilia, Recommendation of the Scientific Subcommittee on Factor VIII and Factor IX of the Scientific and Standardization Committee of the International Society on Thrombosis and Haemostasis, *Thrombosis and Haemostasis*, 2001, 85, p. 560.
8. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
9. "Determination of Coagulation Factor Activities Using the One-Stage Clotting Assay", H48, 2nd edition, CLSI, March 2016.
10. "One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test", Approved Guideline, Second Edition, H47-A2, CLSI, Vol. 28, No. 20, 2008.
11. Andrew, M, Schmidt, B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 991-992, J.B. Lippincott Co., Philadelphia, 1994.
12. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
13. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2012, European Committee for Standardization, 2012.

**SYMBOL LEGEND/ LEGENDE DES SYMBOLES/ LEGENDA DEI
SIMBOLI/ SYMBOL LEGENDE/ LEYENDAS DE SÍMBOLOS /
VERKLARING VAN DE SYMBOLEN 13**

 For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

 Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

 Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

 Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

 Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

 Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenwert / Límite superior de temperatura / Bovenste temperatuurbegrenzing

 Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

 Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

 Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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