VisuDep™-F Frozen Factor X Deficient Plasma

Product #: FRDP100125
FRDP100105

Store at -60°C or lower.

For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuDep-F Frozen Factor X Deficient Plasma is intended for use as a substrate plasma in the quantitative determination of Factor X activity in citrated plasmas using a clot-based assay.

Summary and Principle

Factor X (FX) is a vitamin K-dependent glycoprotein produced in the liver with a plasma concentration of approximately 170 nM. Factor X is expressed as a two-chain molecule with a molecular weight of 59 kDa. The light chain contains a calcium-binding domain consisting of one hydroxy-aspartic acid and 11 γ-carboxyglutamic acid (gla) residues. These residues allow FX to bind to membranes that contain acidic phospholipids in a calcium dependent manner. The heavy chain consists of the catalytic domain, carbohydrate and the activation peptide. Activation of FX to the active enzyme, FXa, results from cleavage at residue Arg⁵² in the heavy chain by a complex of F.IXa, cofactor VIIIa, calcium and negatively charged phospholipid surface (the tenase complex), or by the F.VIIa-tissue factor complex. The FXa generated is a serine protease responsible for the activation of prothrombin to thrombin in the presence of a phospholipid membrane, calcium and cofactor Va. The activity of FXa in plasma is inhibited by antithrombin (ATIII), α₂antitrypsin, α₂macroglobulin and tissue factor pathway inhibitor (TFPI).¹⁻³

Congenital deficiencies of Factor X are rare and are characterized by prolonged prothrombin (PT) and partial thromboplastin times (APTT). Reduced plasma levels of Factor X are observed in patients receiving oral anticoagulant therapy such as warfarin or in patients with vitamin K deficiency or liver disease.⁴

Factor X activity in patient plasma is determined *in vitro* through the use of a modified prothrombin time test (PT). Patient plasma is diluted, mixed with plasma deficient in Factor X and a PT test performed. The result is interpolated using a reference curve generated with dilutions of calibrator plasma mixed with the Factor X deficient plasma. Correction of the clotting time of the deficient plasma is proportional to the activity level of Factor X in the patient plasma.

Reagents

A. Description:

Frozen Factor X immune-depleted plasma with residual Factor X activity of <1%. Remaining coagulation factors have activities of ≥50% and fibrinogen is ≥1.5 g/L.

Box containing 25 x 1 ml vials of frozen Factor X deficient plasma,

Product # FRDP100125

OR

Box containing 5 x 1 ml vials of frozen Factor X deficient plasma,

Product # FRDP100105

B. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this

product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

C. Reagent Preparation/Handling:

Thaw each 1 mL vial for 5 minutes in a 37°C ± 1°C waterbath. Invert gently before use.

D. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma is stable for 8 hours on-board (15-19°C) or at 2-8°C in the original vial. Thawed material should be discarded after use and should not be refrozen.

Specimen Collection and Preparation

Blood is collected into 3.2% sodium Citrate anticoagulant tubes at a ratio of 9 volumes blood to 1 volume anticoagulant and gently mixed by inversion. Refer to CLSI Guideline H21-A5 "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" for further information on specimen collection, handling and storage.⁵

Procedure

A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay reagents e.g. thromboplastin with added calcium chloride, Owren-Koller buffer

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Control Plasmas (e.g. VisuCon-F Frozen Normal and Abnormal Control Plasmas)

Common clinical laboratory equipment and material (pipettes, plastic tubes etc.)

B. Assay Procedure (Manual Method)

Preparation of Calibration Curve Dilutions

i) Prepare calibrator plasma, thromboplastin reagent and assay buffer according to manufacturer's instructions.

ii) Prepare dilutions of the calibrator plasma in plastic tubes using assay buffer to generate a 5-point calibration curve corresponding to 100% to 6.25%, as illustrated in the following table:

Dilution	% Factor X*	Calibrator Plasma Volume	Buffer Volume
1/10	100%	100 µl	900 µl
1/20	50%	500 µl of 1/10 dilution	500 µl
1/40	25%	500 µl of 1/20 dilution	500 µl
1/80	12.5%	500 µl of 1/40 dilution	500 µl
1/160	6.25%	500 µl of 1/80 dilution	500 µl

*Refer to manufacturer's Certificate of Analysis for the lot-specific Factor X activity level of the calibrator plasma.

Preparation of Test Plasma Dilutions

Dilute all test plasma samples 1/10 and 1/20 in plastic tubes using assay buffer. Additional dilutions may be required if the reported results are outside the readable range of the calibration curve or if the test results are not parallel with the calibration curve.⁶

Test Method (Manual)

All calibration plasma and test plasma dilutions should be tested immediately after preparation. Add the following to a test tube or cuvette warmed to 37°C ± 1°C:

- 100 µl of VisuDep-F Factor X deficient plasma
 - 100 µl of calibrator plasma dilution or diluted test sample
- Mix and incubate for 2 minutes at 37°C ± 1°C then add:
- 200 µl of thromboplastin (pre-warmed to 37°C ± 1°C)

Start a timing device the instant the calcium chloride is added to the plasma mixture and record the time to clot formation.

Results

Using log-log graph paper, plot clot times for the calibrator plasma dilutions on the y-axis against the corresponding percent Factor X activity on the x-axis. Draw the best-fit calibration line. Interpolate the Factor X levels of the test samples from this calibration curve and multiply by the appropriate dilution factor. Under the conditions described here, a test sample diluted 1/10 will have a dilution factor of 1 and a dilution of 1/20 will have a dilution factor of 2. Additionally, if the Factor X value for the calibrator plasma is not 100%, the test sample result will need to be multiplied by the appropriate correction factor. For example, if the Factor X level of the calibrator plasma is 105%, multiply the test sample results by 1.05 to obtain the corrected result.

Automatic Test Method

Factor activity assays may be performed using semi-automated or automated coagulation instruments with built in programs for generating the calibration curve and calculating test plasma results. Refer to the coagulation instrument's Operator's Manual for specific assay procedure instructions.

Quality Control

In order to ensure accuracy of the results, normal and abnormal control plasmas such as the VisuCon-F Frozen Normal and Abnormal Controls, should be analyzed regularly. Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Results

Results for test samples and control plasmas can be calculated as described above or using the software program of the coagulation instrument. Ensure values obtained for the control plasmas are within established ranges. If control results are out of range, the Factor X activity values obtained for test samples should be considered suspect.

Limitations of the Procedure

Incorrect sample handling such as insufficient filling of collection tubes, contamination of collection tubes or incorrect type of anticoagulant can lead to erroneous sample results. Use only plasma samples that have been collected into sodium citrate anticoagulant.⁵ Samples that have clotted or are hemolyzed, icteric or lipemic should be avoided. Do not use glass tubes to prepare plasma dilutions as glass may cause activation of the coagulation cascade.⁶⁻⁷

Expected Values

The normal range for Factor X as reported in the literature is typically 70-152%⁸, however, each laboratory should determine a normal range independently. Results measured in 100 healthy individuals indicate a Factor X normal reference range of 80 – 161% (Mean 114%, SD = 20.5).

Performance Characteristics

PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuDep-F Frozen Factor X Deficient Plasma using STA[®] Neoplastine CI+ on an STA Compact. Normal and abnormal control plasmas were tested in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2⁹. The coefficients of variation (% CV) obtained in these precision studies are presented in the table below.

	Within Run	Between Day	Between Run	Within Device
Normal	3.0 – 5.4%	2.7 – 4.4%	0 – 2.6%	5.0 – 7.0%
Abnormal	2.4 – 2.9%	2.4 – 3.4%	2.0 – 3.2%	4.2 – 5.2%

METHOD COMPARISON:

In a method comparison study, 161 test samples containing Factor X ranging from <1% to normal levels were tested using the VisuDep-F Frozen Factor X Deficient Plasma and Precision Biologic's Factor X Deficient Plasma as the reference method. The correlation co-efficient (r) reported in this study was 0.997 (R² = 0.994, y = 1.0086x + 0.2044).

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1. Ichinose A, Davie EW; The Blood Coagulation Factors: Their cDNAs, Genes, and Expression; in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp 19-54, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
2. Steinberg M, Nemerson Y; The Activation of Factor X; in Hemostasis and Thrombosis, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp 91-99, J.B. Lippincott Co., Philadelphia PA, USA, 1982.
3. Ellis V, Scully M, MacGregor I, Kakkar V; Inhibition of Human Factor Xa by Various Plasma Protease Inhibitors; Biochimica et Biophysica Acta 701, pp 24-31, 1982.
4. Roberts, HR, Lefkowitz, JB; Inherited Disorders of Prothrombin Conversion; in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 200 - 218, J.B. Lippincott Co., Philadelphia, PA, USA, 1994.
5. "Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General performance of Coagulation Assays"; Approved Guideline–Fifth Edition. CLSI Document H21-A5, Vol. 28, No. 5. 2008.
6. "Determination of Coagulation Factor Activities Using the One-Stage Clotting Assay", H48, 2nd edition, CLSI, March 2016.
7. "One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test", Approved Guideline, Second Edition, H47-A2, CLSI, Vol. 28, No. 20, 2008.
8. Andrew, M, Schmidt, B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 991-992, J.B. Lippincott Co., Philadelphia, 1994.
9. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
10. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2012, European Committee for Standardization, 2012.

**SYMBOL LEGEND/ LEGENDE DES SYMBOLES/ LEGENDA DEI
SIMBOLI/ SYMBOL LEGENDE/ LEYENDAS DE SÍMBOLOS /
VERKLARING VAN DE SYMBOLEN¹⁰**

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenwert / Límite superior de temperatura / Bovenste temperatuurbegrenzing

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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