****REPRESENTATIVE DATASHEET****
VisuDep™-F Frozen Factor XI Deficient Plasma
Product #: FRDP110125
FRDP110105

Store at -60°C or lower.

 For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuDep-F Frozen Factor XI Deficient Plasma is intended for use as substrate plasma in the quantitative determination of Factor XI activity in citrated plasmas using a clot-based assay.

Summary and Principle

 Factor XI (FXI) is a coagulation protein produced in the liver that circulates in plasma at approximately 5 µg/ml (30 nM). FXI consists of two identical 80 kDa subunits linked by disulphide bonds. Cleavage of FXI by activated factor XII or thrombin converts each subunit into a two-chain form and generates two active sites per FXIa molecule. In plasma, FXI or FXIa circulates in non-covalent 1:1 complex with high molecular weight kininogen, which acts as a cofactor in the activation of FXI by activated factor XII. The activity of FXIa is regulated by platelets and by several proteinase inhibitors including, in order of decreasing importance, C1-inhibitor, α₂antiplasmin, α₁antitrypsin and antithrombin. Heparin has relatively little effect on the rate of inhibition of FXIa by antithrombin. The only known natural substrate for FXIa is factor IX and the only cofactor required for this reaction is ionized calcium¹⁻³.

 Congenital deficiency of Factor XI is a rare disorder that causes mild or moderately severe bleeding, usually after trauma or surgery. Treatment of Factor XI deficiencies can be problematic as clotting activity does not correlate well with clinical severity⁴.

 Factor XI activity in patient plasma is determined *in vitro* through the use of a modified activated partial thromboplastin time test (APTT). Patient plasma is diluted, mixed with plasma deficient in Factor XI and an APTT test performed. The result is interpolated using a reference curve generated with dilutions of calibrator plasma mixed with the Factor XI deficient plasma. Correction of the clotting time of the deficient plasma is proportional to the activity level of Factor XI in the patient plasma.

Reagents
A. Description:

Frozen Factor XI immune-depleted plasma with residual Factor XI activity of <1%. Remaining coagulation factors have activities of ≥50% and fibrinogen is ≥1.5 g/L.

Box containing 25 x 1 ml vials of frozen Factor XI deficient plasma,

Product # FRDP110125

OR

Box containing 5 x 1 ml vials of frozen Factor XI deficient plasma,

Product # FRDP110105
B. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this

product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

C. Reagent Preparation/Handling:

Thaw each 1 mL vial for 5 minutes in a 37°C ± 1°C waterbath. Invert gently before use.

D. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma is stable for 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and not be refrozen.

Specimen Collection and Preparation

 Blood is collected into 3.2% sodium Citrate anticoagulant tubes at a ratio of 9 volumes blood to 1 volume anticoagulant and gently mixed by inversion. Refer to CLSI Guideline H21-A5 "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" for further information on specimen collection, handling and storage.⁵
Procedure
A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay reagents e.g. APTT reagent, calcium chloride, Owren-Koller buffer

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Control Plasmas (e.g. VisuCon-F Frozen Normal and Abnormal Control Plasmas)

Common clinical laboratory equipment and material (pipettes, plastic tubes etc.)

B. Assay Procedure (Manual Method)
Preparation of Calibration Curve Dilutions

i) Prepare calibrator plasma, APTT reagent and assay buffer according to manufacturer's instructions.

ii) Prepare dilutions of the calibrator plasma in plastic tubes using assay buffer to generate a 5-point calibration curve corresponding to 100% to 6.25%, as illustrated in the following table:

Dilution	% Factor XI*	Calibrator Plasma Volume	Buffer Volume
1/10	100%	100 µl	900 µl
1/20	50%	500 µl of 1/10 dilution	500 µl
1/40	25%	500 µl of 1/20 dilution	500 µl
1/80	12.5%	500 µl of 1/40 dilution	500 µl
1/160	6.25%	500 µl of 1/80 dilution	500 µl

*Refer to manufacturer's Certificate of Analysis for the lot-specific Factor XI activity level of the calibrator plasma.

Preparation of Test Plasma Dilutions

 Dilute all test plasma samples 1/10 and 1/20 in plastic tubes using assay buffer. Additional dilutions may be required if the reported results are outside the readable range of the calibration curve or if the test results are not parallel with the calibration curve.⁶
Test Method (Manual)

All calibration plasma and test plasma dilutions should be tested immediately after preparation. Add the following to a test tube or cuvette warmed to 37°C ± 1°C:

- 100 µl of VisuDep-F Factor XI deficient plasma
- 100 µl of calibrator plasma dilution or diluted test sample

Mix and incubate for 1 minute at 37°C ± 1°C then add:

- 100 µl of APTT reagent (pre-warmed to 37°C ± 1°C)

Mix and incubate for 3-5 minutes at 37°C ± 1°C then add:

- 100 µl of calcium chloride (pre-warmed to 37°C ± 1°C)

Start a timing device the instant the calcium chloride is added to the plasma mixture and record the time to clot formation.

Results

Using log-log graph paper, plot clot times for the calibrator plasma dilutions on the y-axis against the corresponding percent Factor XI activity on the x-axis. Draw the best-fit calibration line. Interpolate the Factor XI levels of the test samples from this calibration curve and multiply by the appropriate dilution factor. Under the conditions described here, a test sample diluted 1/10 will have a dilution factor of 1 and a dilution of 1/20 will have a dilution factor of 2. Additionally, if the Factor XI value for the calibrator plasma is not 100%, the test sample result will need to be multiplied by the appropriate correction factor. For example, if the Factor XI level of the calibrator plasma is 105%, multiply the test sample results by 1.05 to obtain the corrected result.

Automatic Test Method

Factor activity assays may be performed using semi-automated or automated coagulation instruments with built in programs for generating the calibration curve and calculating test plasma results. Refer to the coagulation instrument's Operator's Manual for specific assay procedure instructions.

Quality Control

In order to ensure accuracy of the results, normal and abnormal control plasmas such as the VisuCon-F Frozen Normal and Abnormal Controls, should be analyzed regularly. Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Results

Results for test samples and control plasmas can be calculated as described above or using the software program of the coagulation instrument. Ensure values obtained for the control plasmas are within established ranges. If control results are out of range, the Factor XI activity values obtained for test samples should be considered suspect.

Limitations of the Procedure

Incorrect sample handling such as insufficient filling of collection tubes, contamination of collection tubes or incorrect type of anticoagulant can lead to erroneous sample results. Use only plasma samples that have been collected into sodium citrate anticoagulant.⁵ Samples that have clotted or are hemolyzed, icteric or lipemic should be avoided. Do not use glass tubes to prepare plasma dilutions as glass may cause activation of the coagulation cascade.⁶⁻⁷ The presence of lupus anticoagulant in test plasmas can affect the reported factor activity.

Expected Values

The normal range for Factor XI as reported in the literature is typically 67-127%⁸, however, each laboratory should determine a normal range independently. Results measured in 112 healthy individuals indicate a Factor XI normal reference range of 66-145%.

Performance Characteristics

PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuDep-F Frozen Factor XI Deficient Plasma using Hemosil™ APTT-SP reagent on an STA Compact. Normal and abnormal control plasmas were tested in quadruplicate 2 times per day for 20 days as per CLSI guideline EP5-A2⁹. The coefficients of variation (% CV) obtained in these precision studies are presented in the table below.

	Within Run	Between Day	Between Run	Within Device
Normal	5.1 – 6.0%	0 – 5.3%	4.0 – 7.3%	7.0 – 9.3%
Abnormal	5.1 – 5.9%	0.2 – 7.0%	4.5 – 7.2%	7.4 – 10.9%

METHOD COMPARISON:

In a method comparison study 145 test samples containing Factor XI ranging from <1% to normal levels were tested using the VisuDep-F Frozen Factor XI Deficient Plasma and Precision Biologic's Factor XI Deficient Plasma as the reference method. The correlation co-efficient (r) reported in this study was 0.991 ($R^2 = 0.9814$, $y = 1.0137 - 0.1152$).

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Bibliography/ Bibliographie / Bibliografía / Bibliographie / Bibliografía

1. Wuillemin WA, Minnema M, Meijers JCM, Roem D, Erenberg AJM, Nuijens JH, ten Cate H, Hack EC; Inactivation of Factor XIa in Human Plasma Assessed by Measuring Factor XIa-Protease Inhibitor Complexes: Major Role for C1-Inhibitor. *Blood* 85:1517, 1995.
2. DeLa Cadena R, Watchtfoegel YT, Colman RW, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 219-240, J.B. Lippincott Co., Philadelphia, 1994.
3. Baglia FA, Seaman FS, Walsh, PN; The Apple 1 and 4 domains of Factor XI Act to Synergistically Promote the Surface-Mediated Activation of Factor XI by Factor XIIa. *Blood* 85:2078, 1995.
4. Thompson, AR, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 55-80, J.B. Lippincott Co., Philadelphia, 1994.
5. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
6. "Determination of Coagulation Factor Activities Using the One-Stage Clotting Assay", H48, 2nd edition, CLSI, March 2016.
7. "One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test", Approved Guideline, Second Edition, H47-A2, CLSI, Vol. 28, No. 20, 2008.
8. Andrew, M, Schmidt, B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 991-992, J.B. Lippincott Co., Philadelphia, 1994.
9. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
10. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2016, European Committee for Standardization, 2016.

**SYMBOL LEGEND / LEGENDE DES SYMBOLES / LEGENDA DEI
SIMBOLI / SYMBOL LEGENDE / LEYENDAS DE SÍMBOLOS /
VERKLARING VAN DE SYMBOLEN¹⁰**

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenzwert / Límite superior de temperatura / Bovenste temperatuurbegrenzing

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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