VisuDep™-F Frozen Factor XII Deficient Plasma

Product #: FRDP120125
FRDP120105

Store at -60°C or lower.

For *in vitro* diagnostic use.

Intended Use

ENGLISH

The VisuDep-F Frozen Factor XII Deficient Plasma is intended for use as substrate plasma in the quantitative determination of Factor XII activity in citrated plasmas using a clot-based assay.

Summary and Principle

Factor XII (FXII, Hageman factor) is a 76 kDa single chain glycoprotein produced in the liver with a plasma concentration of approximately 400 nM. Upon vascular injury, FXII binds to negatively charged extravascular surfaces, which facilitate activation of the zymogen to the active serine protease. Cleavage of FXII by kallikrein produces the enzyme α FXIIa, consisting of a 28 kDa light chain containing the protease domain, and a 52 kDa heavy chain containing the anionic surface-binding domain. Substrates for surface bound FXIIa include prekallikrein (PK) and factor XI as well as the procofactor high-molecular weight kininogen. The activation of these substrates results in positive feedback activation of FXII. Further cleavage of α FXIIa by kallikrein produces the 28 kDa fragment β FXIIa which has reduced procoagulant activity as it lacks the anionic surface-binding domain, but is capable of fluid-phase activation of PK, factor VII and complement C1. The activity of FXIIa in plasma is regulated predominantly by C1-Inhibitor, with minor contributions by α_2 antiplasmin, α_2 macroglobulin and antithrombin, even in the presence of therapeutic levels of heparin¹⁻³.

Congenital deficiencies of FXII are very rare and are characterized by prolonged activated partial thromboplastin time (APTT) but with no apparent associated bleeding disorder.⁴

Factor XII activity in patient plasma is determined *in vitro* through the use of a modified activated partial thromboplastin time test (APTT). Patient plasma is diluted, mixed with plasma deficient in Factor XII and an APTT test performed. The result is interpolated using a reference curve generated with dilutions of calibrator plasma mixed with the Factor XII deficient plasma. Correction of the clotting time of the deficient plasma is proportional to the activity level of Factor XII in the patient plasma.

Reagents

A. Description:

Frozen Factor XII immune-depleted plasma with residual Factor XII activity of <1%. Remaining coagulation factors have activities of \geq 50% and fibrinogen is \geq 1.5 g/L.

Box containing 25 x 1 ml vials of frozen Factor XII deficient plasma,

Product # FRDP120125

OR

Box containing 5 x 1 ml vials of frozen Factor XII deficient plasma,

Product # FRDP120105

B. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each unit of source plasma used in the preparation of this product has been tested by FDA approved methods and found non-reactive for Hepatitis B surface antigen (HbsAg), syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA. As no test can offer complete assurance that

products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material. The disposal of materials must be carried out according to current local regulations.

C. Reagent Preparation/Handling:

Thaw each 1 mL vial for 5 minutes in a 37°C \pm 1°C waterbath. Invert gently before use.

D. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma is stable for 4 hours on-board (15-19°C) or 8 hours at 2-8°C in original vial. Thawed material should be discarded after use and not be refrozen.

Specimen Collection and Preparation

Blood is collected into 3.2% sodium Citrate anticoagulant tubes at a ratio of 9 volumes blood to 1 volume anticoagulant and gently mixed by inversion. Refer to CLSI Guideline H21-A5 "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" for further information on specimen collection, handling and storage.⁵

Procedure

A. Materials Required (but not provided)

Waterbath

Coagulation Instrument or assay system

Assay reagents e.g. APTT reagent, calcium chloride, Owren-Koller buffer

Calibrator Plasma (e.g. VisuCal-F Frozen Calibrator Plasma)

Control Plasmas (e.g. VisuCon-F Frozen Normal and Abnormal Control Plasmas)

Common clinical laboratory equipment and material (pipettes, plastic tubes etc.)

B. Assay Procedure (Manual Method)

Preparation of Calibration Curve Dilutions

i) Prepare calibrator plasma, APTT reagent and assay buffer according to manufacturer's instructions.

ii) Prepare dilutions of the calibrator plasma in plastic tubes using assay buffer to generate a 5-point calibration curve corresponding to 100% to 6.25%, as illustrated in the following table:

Dilution	% Factor XII*	Calibrator Plasma Volume	Buffer Volume
1/10	100%	100 μ l	900 μ l
1/20	50%	500 μ l of 1/10 dilution	500 μ l
1/40	25%	500 μ l of 1/20 dilution	500 μ l
1/80	12.5%	500 μ l of 1/40 dilution	500 μ l
1/160	6.25%	500 μ l of 1/80 dilution	500 μ l

*Refer to manufacturer's Certificate of Analysis for the lot-specific Factor XII activity level of the calibrator plasma.

Preparation of Test Plasma Dilutions

Dilute all test plasma samples 1/10 and 1/20 in plastic tubes using assay buffer. Additional dilutions may be required if the reported results are outside the readable range of the calibration curve or if the test results are not parallel with the calibration curve.⁶

Test Method (Manual)

All calibration plasma and test plasma dilutions should be tested immediately after preparation. Add the following to a test tube or cuvette warmed to 37°C \pm 1°C:

- 100 μ l of VisuDep-F Factor XII deficient plasma
- 100 μ l of calibrator plasma dilution or diluted test sample

Mix and incubate for 1 minute at 37°C \pm 1°C then add:

- 100 μ l of APTT reagent (pre-warmed to 37°C \pm 1°C)
Mix and incubate for 3-5 minutes at 37°C \pm 1°C then add:
- 100 μ l of calcium chloride (pre-warmed to 37°C \pm 1°C)

Start a timing device the instant the calcium chloride is added to the plasma mixture and record the time to clot formation.

Results

Using log-log graph paper, plot clot times for the calibrator plasma dilutions on the y-axis against the corresponding percent Factor XII activity on the x-axis. Draw the best-fit calibration line. Interpolate the Factor XII levels of the test samples from this calibration curve and multiply by the appropriate dilution factor. Under the conditions described here, a test sample diluted 1/10 will have a dilution factor of 1 and a dilution of 1/20 will have a dilution factor of 2. Additionally, if the Factor XII value for the calibrator plasma is not 100%, the test sample result will need to be multiplied by the appropriate correction factor. For example, if the Factor XII level of the calibrator plasma is 105%, multiply the test sample results by 1.05 to obtain the corrected result.

Automatic Test Method

Factor activity assays may be performed using semi-automated or automated coagulation instruments with built in programs for generating the calibration curve and calculating test plasma results. Refer to the coagulation instrument's Operator's Manual for specific assay procedure instructions.

Quality Control

In order to ensure accuracy of the results, normal and abnormal control plasmas such as the VisuCon-F Frozen Normal and Abnormal Controls, should be analyzed regularly. Each laboratory should establish and maintain its own quality control ranges for each particular instrument-reagent system used. If appropriate control values are not obtained, assess the components of the assay system including reagents, substrate plasmas, calibrator plasma or instrumentation to try and identify potential sources of error.

Results

Results for test samples and control plasmas can be calculated as described above or using the software program of the coagulation instrument. Ensure values obtained for the control plasmas are within established ranges. If control results are out of range, the Factor XII activity values obtained for test samples should be considered suspect.

Limitations of the Procedure

Incorrect sample handling such as insufficient filling of collection tubes, contamination of collection tubes or incorrect type of anticoagulant can lead to erroneous sample results. Use only plasma samples that have been collected into sodium citrate anticoagulant.⁵ Samples that have clotted or are hemolyzed, icteric or lipemic should be avoided. Do not use glass tubes to prepare plasma dilutions as glass may cause activation of the coagulation cascade.⁶⁻⁷ The presence of lupus anticoagulant in test plasmas can affect the reported factor activity.

Expected Values

The normal range for Factor XII as reported in the literature is typically 52-164%⁸; however, each laboratory should determine a normal range independently. Results measured in 101 healthy individuals indicate a Factor XII normal reference range of 53-149%.

Performance Characteristics

PRECISION:

Within-run, between day, between-run and within device precision were assessed for three lots of the VisuDep-F Frozen Factor XII Deficient Plasma using STA[®] CK-Prest on an STA Compact. Normal and abnormal control plasmas were tested in duplicate 2 times per day for 20 days as per CLSI guideline EP5-A2⁹. The coefficients of variation (% CV) obtained in these precision studies are presented in the table below.

	Within Run	Between Day	Between Run	Within Device
Normal	4.5 – 7.1%	0.1 – 3.5%	2.6 – 5.1%	6.2 – 8.6%
Abnormal	5.1 – 6.5%	0.7 – 4.1%	3.2 – 5.3%	7.3 – 7.6%

METHOD COMPARISON:

In a method comparison study 149 test samples containing Factor XII ranging from <1% to normal levels were tested using the VisuDep-F Frozen Factor XII Deficient Plasma and Precision Biologic's Factor XII Deficient Plasma as the reference method. The correlation co-efficient (r) reported in this study was 0.992 ($R^2 = 0.984$, $y = 1.0597x + 0.2023$).

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1. DeLa Cadena R, Watchfogel YT, Colman RW, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 219-240, J.B. Lippincott Co., Philadelphia, 1994.
2. Tankersley DL, Alving BM, Finlayson JS; Preparation of β FXIIa (Hageman Factor Fragment) from Human Plasma. Thrombosis Research 25 pp 307-317, 1982.
3. Pixley RA, Schapira M, Coleman RW; The regulation of Human Factor XIIa by Plasma Proteinase Inhibitors. JBC 260, pp 1723-1729, 1985.
4. Thompson, AR, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 55-80, J.B. Lippincott Co., Philadelphia, 1994.
5. "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline, Fifth Edition. H21-A5, CLSI, Vol. 28. No. 5, 2008.
6. "Determination of Coagulation Factor Activities Using the One-Stage Clotting Assay", H48, 2nd edition, CLSI, March 2016.
7. "One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test", Approved Guideline, Second Edition, H47-A2, CLSI, Vol. 28, No. 20, 2008.
8. Andrew, M, Schmidt, B, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 991-992, J.B. Lippincott Co., Philadelphia, 1994.
9. "Evaluation of Precision Performance of Quantitative Measurement Methods" Approved Guideline, Second Edition, EP5-A2, CLSI, Vol. 24, No. 25, 2004.
10. "Medical Devices. Symbols to be used with medical device labels, labeling and information to be supplied. General Requirements", EN ISO 15223-1:2012, European Committee for Standardization, 2012.

**SYMBOL LEGEND / LEGENDE DES SYMBOLES / LEGENDA DEI
SIMBOLI / SYMBOL LEGENDE / LEYENDAS DE SÍMBOLOS /
VERKLARING VAN DE SYMBOLEN¹⁰**

For in vitro diagnostic use / Pour utilisation diagnostique in vitro / Per uso diagnostico in vitro / für In-Vitro-Diagnose / Para uso en diagnóstico in vitro / Voor in-vitrodiagnostiek

Batch code / Code du lot / Codice lotto / Chargencode / Código de lote / Partnummer

Authorized Representative / Représentant agréé / Rappresentante autorizzato / Autorisierter Vertreter / Representante autorizado / Gemachtigde Vertegenwoordiger

Expiry date / Date d'expiration / Data di scadenza / Verfallsdatum / Fecha de caducidad / Vervaldatum

Catalogue Number / Référence du catalogue / Numero di catalogo/ Katalognummer / Número de catálogo / Catalogusnummer

Upper limit of temperature / Limite supérieure de température / Limite superiore di temperatura / Oberer Temperaturgrenzwert / Límite superior de temperatura / Bovenste temperatuurbegrenzing

Manufacturer/ Fabricant / Produttore / Hersteller / Fabricante / Fabrikant

Consult instructions for use / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Gebrauchsanleitung konsultieren / Consultar las instrucciones de uso / Raadpleeg de gebruiksaanwijzing

Biological Risks / Risques biologiques / Rischio Biologico / Biologische Risiken / Riesgos biológicos / Biologisch risico

Limited Warranty:

This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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