



Product Data Sheet for:  
**Non-Immune Rabbit Immunoglobulin  
Purified IgG**

**Catalogue #:** NIR-IG  
**Lot #:** XXXX  
**Expiry:** XXXX

**Presentation:** XXXX ml of purified immunoglobulin from rabbit serum. The IgG concentration is XXXX mg/ml in 10 mM HEPES, pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v). Total protein is **10 mg**.

**Storage:** Store between -10 and -20°C. Product will become viscous but will not freeze. Avoid storage in frost-free freezers. Keep vial tightly capped. Allow product to warm to room temperature and gently mix before use.

**Specificity:** No reaction was observed with normal human plasma using immunoelectrophoresis and immunodiffusion techniques.

**Application:** This product is suitable as a negative control in procedures where rabbit antibodies are used.

***For Research Use Only***

Related products:

Product: <b>NIS-IG</b>	Non-immune sheep immunoglobulin, purified IgG	20 mg vial
Product: <b>NIS-FITC</b>	Non-immune sheep immunoglobulin, FITC-conjugated purified IgG	0.1 mg vial
Product: <b>NIS-BIO</b>	Non-immune sheep immunoglobulin, biotinylated purified IgG	0.1 mg vial
Product: <b>NIG-IG</b>	Non-immune goat immunoglobulin, purified IgG	10 mg vial
Product: <b>NIG-FITC</b>	Non-immune goat immunoglobulin, FITC-conjugated purified IgG	0.1 mg vial
Product: <b>NIG-BIO</b>	Non-immune goat immunoglobulin, biotinylated purified IgG	0.1 mg vial

Visit our website ([www.affinitybiologicals.com](http://www.affinitybiologicals.com)) for details.

**Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory**

**Limited Warranty:** This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.