



Product Data Sheet for:  
**Non-Immune Sheep Immunoglobulin  
FITC-Conjugated Purified IgG**

**Catalogue #:** NIS-FITC  
**Lot #:** XXXX  
**Expiry:** XXXX

**Presentation:** XXXX ml of purified immunoglobulin from sheep serum conjugated to fluorescein isothiocyanate (FITC). The IgG-FITC concentration is XXXX mg/ml in phosphate-buffered saline (PBS) containing 1 mg/ml bovine albumin and 0.1% sodium azide (w/v), pH 7.4. The incorporation of FITC is XXXX moles fluorescein per mole IgG as determined spectrophotometrically. Total protein is 0.1 mg.

**Storage:** Store at 2-8°C and protect from light.

**Specificity:** Prior to conjugation, no reaction was observed with normal human plasma using immunoelectrophoresis and immunodiffusion techniques.

**Application:** This product is suitable as a negative control in procedures where sheep antibodies are used.

***For Research Use Only***

Related products:

Product: <b>NIS-IG</b>	Non-immune sheep immunoglobulin, purified IgG	20 mg vial
Product: <b>NIS-BIO</b>	Non-immune sheep immunoglobulin, biotinylated purified IgG	0.1 mg vial
Product: <b>NIG-IG</b>	Non-immune goat immunoglobulin, purified IgG	10 mg vial
Product: <b>NIG-FITC</b>	Non-immune goat immunoglobulin, FITC-conjugated purified IgG	0.1 mg vial
Product: <b>NIG-BIO</b>	Non-immune goat immunoglobulin, biotinylated purified IgG	0.1 mg vial
Product: <b>NIR-IG</b>	Non-immune rabbit immunoglobulin, purified IgG	10 mg vial

Visit our website ([www.affinitybiologicals.com](http://www.affinitybiologicals.com)) for details.

**Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory**

**Limited Warranty:** This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.