



Product Data Sheet for:
Donkey anti-sheep IgG (Heavy & Light Chain)
Biotinylated Affinity-Purified IgG

Catalogue #: DAS-APBIO

Lot #: XXXX

Expiry: XXXX

Immunogen: Whole IgG (heavy and light chain) purified from sheep serum.

Presentation: XXXX ml of donkey immunoglobulin affinity purified on immobilized sheep IgG and conjugated to biotin. The APIgG-biotin concentration is XXXX mg/ml in phosphate-buffered saline (PBS) containing 1 mg/ml bovine albumin and 0.1% sodium azide (w/v), pH 7.4. The incorporation of biotin is XXXX moles biotin per mole IgG as determined by HABA assay. Total protein is 0.5 mg.

Storage: Store at 2-8°C.

Specificity: Immunoelectrophoresis of this antibody prior to conjugation gave a single precipitin line of identical mobility to that obtained with purified sheep IgG. There is 100% cross-reactivity with goat IgG when tested by ELISA.

Application: This product is suitable as a source of biotinylated antibodies to sheep IgG.

For Research Use Only

Related products:

Product: GAM-APHRP	Goat anti-Mouse, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: SAR-APHRP	Sheep anti-Rabbit, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: DAS-APHRP	Donkey anti-sheep, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: DAS-APFTC	Donkey anti-sheep, FITC-conjugated affinity purified IgG	0.5 mg vial

Visit our website (www.affinitybiologicals.com) for details.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Limited Warranty: This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.