



Product Data Sheet for:  
**Donkey anti-sheep IgG (Heavy & Light Chain)**  
**FITC-Conjugated Affinity-Purified IgG**

**Catalogue #:** DAS-APFTC

**Lot #:** XXXX

**Expiry:** XXXX

**Immunogen:** Whole IgG (heavy and light chain) purified from sheep serum.

**Presentation:** XXXX ml of donkey immunoglobulin affinity purified on immobilized sheep IgG and conjugated to fluorescein isothiocyanate (FITC). The APIgG-FITC concentration is XXXX mg/ml in phosphate-buffered saline (PBS) containing 1 mg/ml bovine albumin and 0.1% sodium azide (w/v), pH 7.4. The incorporation of FITC is XXXX moles fluorescein per mole IgG as determined spectrophotometrically. Total protein is 0.5 mg.

**Storage:** Store at 2-8°C and protect from light.

**Specificity:** Immuno-electrophoresis of this antibody prior to conjugation gave a single precipitin line of identical mobility to that obtained with purified sheep IgG. There is 100% cross-reactivity with goat IgG when tested by ELISA.

**Application:** This product is suitable as a source of FITC-conjugated antibodies to sheep IgG.

***For Research Use Only***

Related products:

Product: <b>GAM-APHRP</b>	Goat anti-Mouse, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: <b>SAR-APHRP</b>	Sheep anti-Rabbit, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: <b>DAS-APHRP</b>	Donkey anti-sheep, peroxidase conjugated affinity-purified IgG	0.5 mg vial
Product: <b>DAS-APBIO</b>	Donkey anti-sheep, biotinylated affinity purified IgG	0.5 mg vial

Visit our website ([www.affinitybiologicals.com](http://www.affinitybiologicals.com)) for details.

**Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory**

**Limited Warranty:** This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.