



Affinity Biologicals™ INC.

1348 Sandhill Drive, Ancaster, Ontario, Canada L9G 4V5
905-304-9896 • 800-903-6020 • fax 905-304-9897

****REPRESENTATIVE DATASHEET****

VisuCon™ Normal Donor Set

Product #: NDSET

Store at -60°C or lower.

For Research Use Only. Not for use in diagnostic procedures.

Product Description

The VisuCon™ Normal Donor Set is a set of individual normal donor plasmas intended for use in establishing normal ranges for assays where required. This product is intended for Research Use Only and must not be used for diagnostic applications.

The plasmas of the VisuCon™ Normal Donor Set are citrated, platelet poor plasmas collected from male and female donors between the ages of 18 and 65 years who are non-smokers and medication-free. All plasmas are screened and found negative for FV Leiden (APCR assay) and Lupus Anticoagulant and report results within the normal range for Prothrombin time (PT), Activated Partial Thromboplastin Time (APTT), fibrinogen (clottable) and Factor II activity.

Materials Provided

A. Description:

Box containing 1 mL vial of plasma from each of 25 screened normal human donors. Custom size packaging also available for this product.

The VisuCon™ Normal Donor Set contains a lot-specific Certificate of Analysis indicating the age and gender of each normal donor plasma.

B. Reagent Preparation/Handling:

Thaw each vial at 37°C (+/- 1°C) in a waterbath for 5 minutes. Invert gently before use.

C. Storage and Stability:

Vials are stable until the expiration date stated on the vial when stored at -60°C or lower. Once thawed, plasma should be used immediately. Discard thawed plasma after use.

D. Precautions and Warnings:

This product is intended for use by personnel trained in laboratory procedures and universal precautions for the use of potentially biohazardous substances. This product contains human source material. Each donor plasma used in the preparation of this product has been tested by FDA approved methods and found negative for the presence of Human Immunodeficiency Virus (HIV-1/2, rDNA), Hepatitis B surface antigen (HbsAg) as well as for Hepatitis C (HCV) and STS. As no test can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious material.

Limited Warranty: This product is warranted to perform in accordance with its labeling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.

Manufactured in Canada by:
AFFINITY BIOLOGICALS INC.
1348 Sandhill Drive
Ancaster ON CANADA L9G 4V5
Tel: (905) 304-9896
(800) 903-6020
Fax: (905) 304-9897
info@affinitybiologicals.com